

KD-Validated Anti-TNF Alpha Induced Protein 8 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1726

Product Overview


Name	KD-Validated Anti-TNF Alpha Induced Protein 8 Rabbit Monoclonal Antibody
Accession(Primary)	O95379
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TNFAIP8
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 23 kDa , observed , 18,21 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (O95379)


Gene Name	TNFAIP8
About Gene Symbol	TNFAIP8 (TNF Alpha Induced Protein 8) is a cytokine and immune signaling protein involved in cell migration and invasion. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC