

KD-Validated Anti-SLC25A12 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1727

Product Overview


Name	KD-Validated Anti-SLC25A12 Rabbit Monoclonal Antibody
Accession(Primary)	O75746
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SLC25A12
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 75 kDa , observed , 72 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (O75746)


Gene Name	SLC25A12
About Gene Symbol	SLC25A12 (Solute Carrier Family 25 Member 12) is a membrane transporter involved in membrane transport. It is an important research target in pharmacology, drug resistance, and membrane biology and is widely studied to elucidate cellular mechanisms and id



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC