

KD-Validated Anti-Somatostatin Receptor 5 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1734

Product Overview

Name	KD-Validated Anti-Somatostatin Receptor 5 Rabbit Monoclonal Antibody
Accession(Primary)	P35346
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SSTR5
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 39 kDa , observed , 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000


Target information (P35346)

Gene Name	SSTR5
About Gene Symbol	SSTR5 (Somatostatin Receptor 5) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-somatostatin receptor 5 antibody (Cat#ABG1734). Total cell lysates (5 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-somatostatin receptor 5 antibody (Cat#ABG1734, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




©2024 Genuin Biotechnol

Western blotting analysis using anti-Somatostatin receptor 5 antibody (Cat#ABG1734). Somatostatin receptor 5 expression in wild type (WT) and Somatostatin receptor 5 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Somatostatin receptor 5 antibody (Cat#ABG1734, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Somatostatin receptor 5 expression in HepG2 cells using anti-Somatostatin receptor 5 antibody (Cat#ABG1734, 1:2,000). Green, isotype control; red, Somatostatin receptor 5.



Immunocytochemical staining of HepG2 cells with anti-Somatostatin receptor 5 antibody (Cat#ABG1734, 1:1,000). Nuclei were stained blue with DAPI; Somatostatin receptor 5 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.