

KD-Validated Anti-TIMM50 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1735

Product Overview

Name	KD-Validated Anti-TIMM50 Rabbit Monoclonal Antibody
Accession(Primary)	Q3ZCQ8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TIMM50
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 40 kDa , observed , 36 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q3ZCQ8)

Gene Name	TIMM50
About Gene Symbol	TIMM50 (Translocase Of Inner Mitochondrial Membrane 50) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



