

KD-Validated Anti-TIA1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1738

Product Overview

Name	KD-Validated Anti-TIA1 Rabbit Monoclonal Antibody
Accession(Primary)	P31483
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Species Reactivity	Human, Mouse, Rat
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Target information (P31483)

Gene Name	TIA1
About Gene Symbol	TIA1 (TIA1 Cytotoxic Granule Associated RNA Binding Protein) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to



