

KD-Validated Anti-Sorting Nexin 5 Rabbit Monoclonal

Angiody

Rabbit monoclonal antibody

ABG1740

Product Overview

Name	KD-Validated Anti-Sorting Nexin 5 Rabbit Monoclonal Antibody
Accession(Primary)	Q9Y5X3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SNX5
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 47 kDa , observed , 47 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q9Y5X3)

Gene Name	SNX5
About Gene Symbol	SNX5 (Sorting Nexin 5) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



