

KD-Validated Anti-Transmembrane Protein 43 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1741

Product Overview


Name	KD-Validated Anti-Transmembrane Protein 43 Rabbit Monoclonal Antibody
Accession(Primary)	Q9BTV4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TMEM43
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 45 kDa , observed , 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9BTV4)


Gene Name	TMEM43
About Gene Symbol	TMEM43 (Transmembrane Protein 43) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



©2024 Genuin Biotechnol



Copyright ©2024 Genuin Biotechnologies LLC