

KD-Validated Anti-Ferritin heavy chain 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1742

Product Overview

Name	KD-Validated Anti-Ferritin heavy chain 1 Rabbit Monoclonal Antibody
Accession(Primary)	P02794
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Species Reactivity	Human
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Target information (P02794)

Gene Name	FTH1
About Gene Symbol	FTH1 (Ferritin Heavy Chain 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

