

KD-Validated Anti-VPS26 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1745


Product Overview

Name	KD-Validated Anti-VPS26 Rabbit Monoclonal Antibody
Accession(Primary)	O75436
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human VPS26
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 38 kDa , observed , 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (O75436)

Gene Name	VPS26A
About Gene Symbol	VPS26A (VPS26 Retromer Complex Component A) is a protein kinase and signaling regulator involved in vesicular trafficking. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-VPS26 antibody (Cat#ABG1745, 1:1,000).

Nuclei were stained blue with DAPI; VPS26 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.