

KD-Validated Anti-TNF Receptor Associated Protein 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1749

Product Overview


Name	KD-Validated Anti-TNF Receptor Associated Protein 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q12931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TRAP1
Species Reactivity	Human
Calculated MW	Predicted, 80 kDa , observed , 75 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q12931)


Gene Name	TRAP1
About Gene Symbol	TRAP1 (TNF Receptor Associated Protein 1) is a molecular chaperone involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-TNF receptor associated protein 1 antibody (Cat#ABG1749, 1:1,000). Nuclei were stained blue with DAPI; TNF receptor associated protein 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.