

KD-Validated Anti-Synaptosome Associated Protein 25

Rabbit monoclonal antibody


ABG1757

Product Overview


Name	KD-Validated Anti-Synaptosome Associated Protein 25 Rabbit Monoclonal Antibody
Accession(Primary)	P60880
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SNAP25
Species Reactivity	Human
Calculated MW	Predicted, 23 kDa , observed , 25 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000

Target information (P60880)

Gene Name	SNAP25
About Gene Symbol	SNAP25 (Synaptosome Associated Protein 25) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



2024 Genuin Biotechn