

KD-Validated Anti-ZPR1 Zinc Finger Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1758

Product Overview


Name	KD-Validated Anti-ZPR1 Zinc Finger Rabbit Monoclonal Antibody
Accession(Primary)	O75312
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ZNF259
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 51 kDa , observed , 51 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (O75312)


Gene Name	ZPR1
About Gene Symbol	ZPR1 (ZPR1 Zinc Finger) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-ZPR1 zinc finger antibody (Cat#ABG1758, 1:1,000). Nuclei were stained blue with DAPI; ZPR1 zinc finger was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 µm.