

KD-Validated Anti-UBE2L3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1759


Product Overview

Name	KD-Validated Anti-UBE2L3 Rabbit Monoclonal Antibody
Accession(Primary)	P68036
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human UBE2L3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 18 kDa , observed , 17 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P68036)

Gene Name	UBE2L3
About Gene Symbol	UBE2L3 (Ubiquitin Conjugating Enzyme E2 L3) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-UBE2L3 antibody (Cat#ABG1759, 1:1,000).

Nuclei were stained blue with DAPI; UBE2L3 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.