

KD-Validated Anti-TATA Element Modulatory Factor 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1762


Product Overview

Name	KD-Validated Anti-TATA Element Modulatory Factor 1 Rabbit Monoclonal Antibody
Accession(Primary)	P82094
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TMF
Species Reactivity	Human
Calculated MW	Predicted, 123 kDa , observed , 150 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry(FCM):1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P82094)

Gene Name	TMF1
About Gene Symbol	TMF1 (TATA Element Modulatory Factor 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-TMF antibody (Cat#ABG1762, 1:1,000). Nuclei were stained blue with DAPI; TMF was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.