

KD-Validated Anti-THO Complex Subunit 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1763

Product Overview


Name	KD-Validated Anti-THO Complex Subunit 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q96FV9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Nuclear Matrix Protein p84
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 76 kDa , observed , 76 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q96FV9)


Gene Name	THOC1
About Gene Symbol	THOC1 (THO Complex Subunit 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



right ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-THO complex subunit 1 antibody (Cat#ABG1763, 1:1,000). Nuclei were stained blue with DAPI; THO complex subunit 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.