

KD-Validated Anti-Sorting Nexin 4 Rabbit Monoclonal

Angiody

Rabbit monoclonal antibody

ABG1764

Product Overview

Name	KD-Validated Anti-Sorting Nexin 4 Rabbit Monoclonal Antibody
Accession(Primary)	O95219
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SNX4
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 52 kDa , observed , 52 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (O95219)

Gene Name	SNX4
About Gene Symbol	SNX4 (Sorting Nexin 4) is an autophagy regulatory protein involved in autophagy. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

