

KD-Validated Anti-Integrin Subunit Alpha 5 Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG1770

Product Overview


Name	KD-Validated Anti-Integrin Subunit Alpha 5 Rabbit Polyclonal Antibody
Accession(Primary)	P08648
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CD49e LC
Species Reactivity	Human
Calculated MW	Predicted, 115 kDa , observed , 18 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000

Target information (P08648)

Gene Name	ITGA5
About Gene Symbol	ITGA5 (Integrin Subunit Alpha 5) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC