

KD-Validated Anti-Integrin Subunit Alpha V Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG1771

Product Overview


Name	KD-Validated Anti-Integrin Subunit Alpha V Rabbit Polyclonal Antibody
Accession(Primary)	P06756
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CD51 HC
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 116 kDa , observed , 125 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000

Target information (P06756)

Gene Name	ITGAV
About Gene Symbol	ITGAV (Integrin Subunit Alpha V) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC