

KD-Validated Anti-Proliferating Cell Nuclear Antigen Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG1774

Product Overview


Name	KD-Validated Anti-Proliferating Cell Nuclear Antigen Rabbit Polyclonal Antibody
Accession(Primary)	P12004
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PCNA
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 29 kDa , observed , 36 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P12004)

Gene Name	PCNA
About Gene Symbol	PCNA (Proliferating Cell Nuclear Antigen) is a protein kinase and signaling regulator involved in cell cycle regulation and adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC