

KD-Validated Anti-EEF2K Rabbit Polyclonal Antibody

Rabbit polyclonal antibody

ABG1776

Product Overview

Name	KD-Validated Anti-EEF2K Rabbit Polyclonal Antibody
Accession(Primary)	O00418
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human EEF2K
Species Reactivity	Human, Rat
Calculated MW	Predicted, 82 kDa , observed , 105 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000

Target information (O00418)

Gene Name	EEF2K
About Gene Symbol	EEF2K (Eukaryotic Elongation Factor 2 Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to

