

KD-Validated Anti-NR3C1 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1779

Product Overview

Name	KD-Validated Anti-NR3C1 Mouse Monoclonal Antibody
Accession(Primary)	P04150
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant protein of human Glucocorticoid Receptor
Species Reactivity	Human
Calculated MW	Predicted, 86 kDa , observed , 94 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000


Target information (P04150)

Gene Name	NR3C1
About Gene Symbol	NR3C1 (Nuclear Receptor Subfamily 3 Group C Member 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Glucocorticoid Receptor antibody (Cat#ABG1779). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Glucocorticoid Receptor antibody (Cat#ABG1779, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Glucocorticoid Receptor antibody (Cat#ABG1779). Glucocorticoid Receptor expression in wild type (WT) and Glucocorticoid Receptor shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Glucocorticoid Receptor antibody (Cat#ABG1779, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.