

KD-Validated Anti-LAMP1 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1781

Product Overview

Name	KD-Validated Anti-LAMP1 Mouse Monoclonal Antibody
Accession(Primary)	P11279
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant protein of human CD107a
Species Reactivity	Human
Calculated MW	Predicted, 45 kDa , observed , 90-120 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000

Target information (P11279)

Gene Name	LAMP1
About Gene Symbol	LAMP1 (Lysosomal Associated Membrane Protein 1) is a cell surface marker involved in adaptive immunity and vesicular trafficking. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

