

KD-Validated Anti-Cyclin Dependent Kinase 4 Mouse Monoclonal Antibody


ABG1782

Product Overview


Name	KD-Validated Anti-Cyclin Dependent Kinase 4 Mouse Monoclonal Antibody
Accession(Primary)	P11802
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a kappa
Immunogen	Recombinant protein of human CDK4
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 34 kDa , observed , 33 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P11802)


Gene Name	CDK4
About Gene Symbol	CDK4 (Cyclin Dependent Kinase 4) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC