

KD-Validated Anti-PCNA Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1784

Product Overview

Name	KD-Validated Anti-PCNA Mouse Monoclonal Antibody
Accession(Primary)	P12004
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	A synthesized peptide derived from human PCNA
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 29 kDa , observed , 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P12004)

Gene Name	PCNA
About Gene Symbol	PCNA (Proliferating Cell Nuclear Antigen) is a protein kinase and signaling regulator involved in cell cycle regulation and adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

