

KD-Validated Anti-Exo/Endonuclease G Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1786

Product Overview


Name	KD-Validated Anti-Exo/Endonuclease G Rabbit Monoclonal Antibody
Accession(Primary)	Q9Y2C4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ENDOGL1 / ENGL
Species Reactivity	Human
Calculated MW	Predicted, 4 1 kDa ; Observed, 4 1 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q9Y2C4)


Gene Name	EXOG
About Gene Symbol	EXOG (Exo/Endonuclease G) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Gnuin Biotechnologies LLC



©2024 Gnuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with Exo/Endonuclease G antibody (Cat#ABG1786, 1:1,000). Nuclei were stained blue with DAPI; Exo/Endonuclease G was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 µm.