

KD-Validated Anti-Phosphoglycerate Kinase 1 Mouse Monoclonal Antibody

ABG1789

Product Overview

Name	KD-Validated Anti-Phosphoglycerate Kinase 1 Mouse Monoclonal Antibody
Accession(Primary)	P00558
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG
Immunogen	Recombinant protein of human PGK1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 45 kDa, observed, 44 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P00558)

Gene Name	PGK1
About Gene Symbol	PGK1 (Phosphoglycerate Kinase 1) is a protein kinase and signaling regulator involved in cell migration and invasion. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

