

KD-Validated Anti-Adenylate Kinase 1 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1790

Product Overview



Name	KD-Validated Anti-Adenylate Kinase 1 Mouse Monoclonal Antibody
Accession(Primary)	P00568
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant protein of human Adenylate Kinase 1
Species Reactivity	Human
Calculated MW	Predicted, 45 kDa, observed, 44 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P00568)


Gene Name	AK1
About Gene Symbol	AK1 (Adenylate Kinase 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC