

KD-Validated Anti-Catalase Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1803

Product Overview

Name	KD-Validated Anti-Catalase Rabbit Monoclonal Antibody
Accession(Primary)	P04040
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Catalase
Species Reactivity	Human
Calculated MW	Predicted, 60 kDa, observed, 60 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P04040)

Gene Name	CAT
About Gene Symbol	CAT (Catalase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

