

KD-Validated Anti-Phospho-Cyclin B1 (Ser126) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1806

Product Overview

Name	KD-Validated Anti-Phospho-Cyclin B1 (Ser126) Rabbit Monoclonal Antibody
Accession(Primary)	P14635
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-Cyclin B1 (Ser126)
Species Reactivity	Human
Calculated MW	Predicted, 48 kDa, observed, 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P14635)

Gene Name	CCNB1
About Gene Symbol	CCNB1 (Cyclin B1) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

