

KD-Validated Anti-Growth associated protein 43 Rabbit

Rabbit monoclonal antibody


ABG1808

Product Overview


Name	KD-Validated Anti-Growth associated protein 43 Rabbit Monoclonal Antibody
Accession(Primary)	P17677
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GAP43
Species Reactivity	Human
Calculated MW	Predicted, 25 kDa , observed, 4 8 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:50,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P17677)


Gene Name	GAP43
About Gene Symbol	GAP43 (Growth Associated Protein 43) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HAP-1 cells with Growth associated protein 43 antibody (Cat#ABG1808, 1:1,000). Nuclei were stained blue with DAPI; Growth associated protein 43 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.