

KD-Validated Anti-GATA binding protein 3 Rabbit Monoclonal Antibody


ABG1810

Product Overview


Name	KD-Validated Anti-GATA binding protein 3 Rabbit Monoclonal Antibody
Accession(Primary)	P23771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GATA3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 48 kDa, observed, 40 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P23771)

Gene Name	GATA3
About Gene Symbol	GATA3 (GATA Binding Protein 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genin Biotechnologies LLC



©2024 Genin Biotechnologies LLC