

KD-Validated Anti-Fibroblast Growth Factor 2 Rabbit Monoclonal Antibody



ABG1817


Product Overview

Name	KD-Validated Anti-Fibroblast Growth Factor 2 Rabbit Monoclonal Antibody
Accession(Primary)	P09038
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human FGF2
Species Reactivity	Human
Calculated MW	Predicted, 31 kDa, observed, 18-24kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P09038)

Gene Name	FGF2
About Gene Symbol	FGF2 (Fibroblast Growth Factor 2) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Flow cytometric analysis of FGF2 expression in HepG2 cells using anti-FGF2 antibody (Cat#ABG1817, 1:2,000). Green, isotype control; red, FGF2.