

# KD-Validated Anti-RAB5A Rabbit Monoclonal Antibody

*Rabbit monoclonal antibody*

ABG1828

## Product Overview

|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | KD-Validated Anti-RAB5A Rabbit Monoclonal Antibody                        |
| <b>Accession(Primary)</b>   | P20339                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | Rabbit IgG                                                                |
| <b>Immunogen</b>            | A synthesized peptide derived from human Rab5                             |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Predicted, 24 kDa, observed, 24 kDa                                       |
| <b>Storage Buffer</b>       | Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide. |
| <b>Recommended Dilution</b> | Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000      |

## Target information (P20339)

|                          |                                                                                                                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | RAB5A                                                                                                                                                                                                                                                    |
| <b>About Gene Symbol</b> | RAB5A (RAB5A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |

