

# KD-Validated Anti-BCL10 Immune Signaling Adaptor Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1832

## Product Overview


|                             |                                                                                                              |
|-----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                 | KD-Validated Anti-BCL10 Immune Signaling Adaptor Rabbit Monoclonal Antibody                                  |
| <b>Accession(Primary)</b>   | O95999                                                                                                       |
| <b>Host Species</b>         | Rabbit                                                                                                       |
| <b>Clonality</b>            | Monoclonal                                                                                                   |
| <b>Isotype</b>              | Rabbit IgG                                                                                                   |
| <b>Immunogen</b>            | A synthesized peptide derived from human Bcl10                                                               |
| <b>Species Reactivity</b>   | Human                                                                                                        |
| <b>Calculated MW</b>        | Predicted, 26 kDa, observed, 26 kDa                                                                          |
| <b>Storage Buffer</b>       | Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.                                    |
| <b>Recommended Dilution</b> | Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000 |

## Target information (O95999)


|                          |                                                                                                                                                                                                                                                                            |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | BCL10                                                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | BCL10 (BCL10 Immune Signaling Adaptor) is a cysteine protease (caspase) involved in NF- $\kappa$ B signaling and apoptosis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |




Western blotting analysis using anti-BCL10 antibody (Cat#ABG1832). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-BCL10 antibody (Cat#ABG1832, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-BCL10 antibody (Cat#ABG1832). BCL10 expression in wild type (WT) and BCL10 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-BCL10 antibody (Cat#ABG1832, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of BCL10 expression in HepG2 cells using anti-BCL10 antibody (Cat#ABG1832, 1:2,000). Green, isotype control; red, BCL10.



Immunocytochemical staining of HepG2 cells with anti-BCL10 antibody (Cat#ABG1832, 1:1,000).

Nuclei were stained blue with DAPI; BCL10 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20  $\mu$ m.