

KD-Validated Anti-STAT3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1841

Product Overview


Name	KD-Validated Anti-STAT3 Rabbit Monoclonal Antibody
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human STAT3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 88 kDa, observed, 79 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P40763)


Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-STAT3 antibody (Cat#ABG1841). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-STAT3 antibody (Cat#ABG1841, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-STAT3 antibody (Cat#ABG1841). STAT3 expression in wild type (WT) and STAT3 shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-STAT3 antibody (Cat#ABG1841, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of STAT3 expression in HepG2 cells using anti-STAT3 antibody (Cat#ABG1841, 1:2,000). Green, isotype control; red, STAT3.