

KD-Validated Anti-TANK Binding Kinase 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1844

Product Overview

Name	KD-Validated Anti-TANK Binding Kinase 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q9UHD2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human TBK1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 84 kDa, observed, 74 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9UHD2)

Gene Name	TBK1
About Gene Symbol	TBK1 (TANK Binding Kinase 1) is a serine/threonine protein kinase involved in NF- κ B signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

