

KD-Validated Anti-Cyclin E2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1855

Product Overview

Name	KD-Validated Anti-Cyclin E2 Rabbit Monoclonal Antibody
Accession(Primary)	O96020
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cyclin E2
Species Reactivity	Human
Calculated MW	Predicted, 47 kDa, observed, 47 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (O96020)

Gene Name	CCNE2
About Gene Symbol	CCNE2 (Cyclin E2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

