

KD-Validated Anti-Cofilin 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1858

Product Overview


Name	KD-Validated Anti-Cofilin 1 Rabbit Monoclonal Antibody
Accession(Primary)	P23528
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cofilin
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 19 kDa,observed, 19 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P23528)


Gene Name	CFL1
About Gene Symbol	CFL1 (Cofilin 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC