

KD-Validated Anti-LDL Receptor Related Protein 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1869

Product Overview

Name	KD-Validated Anti-LDL Receptor Related Protein 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q07954
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human LRP1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 505 kDa, observed, 85 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q07954)

Gene Name	LRP1
About Gene Symbol	LRP1 (LDL Receptor Related Protein 1) is a formin actin nucleator involved in TGF- β signaling and adaptive immunity. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

