

KD-Validated Anti-Mitotic Arrest Deficient 2 Like 2 Rabbit Monoclonal Antibody

ABG1870

Product Overview

Name	KD-Validated Anti-Mitotic Arrest Deficient 2 Like 2 Rabbit Monoclonal Antibody
Accession(Primary)	Q9UI95
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Mad2L2
Species Reactivity	Human
Calculated MW	Predicted, 24 kDa, observed, 24 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9UI95)

Gene Name	MAD2L2
About Gene Symbol	MAD2L2 (Mitotic Arrest Deficient 2 Like 2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

