

KD-Validated Anti-PCYT1A Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1875

Product Overview

Name	KD-Validated Anti-PCYT1A Rabbit Monoclonal Antibody
Accession(Primary)	P49585
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CTA
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 42 kDa, observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P49585)

Gene Name	PCYT1A
About Gene Symbol	PCYT1A (Phosphate Cytidyltransferase 1A, Choline) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanism

