

KD-Validated Anti-PRMT1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1876

Product Overview

Name	KD-Validated Anti-PRMT1 Rabbit Monoclonal Antibody
Accession(Primary)	Q99873
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PRMT1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 42 kDa,observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (Q99873)

Gene Name	PRMT1
About Gene Symbol	PRMT1 (Protein Arginine Methyltransferase 1) is a methyltransferase involved in JAK-STAT/cytokine signaling and chromatin remodeling and epigenetics. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

