

KD-Validated Anti-RELB Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1878

Product Overview

Name	KD-Validated Anti-RELB Rabbit Monoclonal Antibody
Accession(Primary)	Q01201
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Rel B
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 62 kDa, observed, 70 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q01201)

Gene Name	RELB
About Gene Symbol	RELB (RELB Proto-Oncogene, NF- κ B Subunit) is a proto-oncogene involved in NF- κ B signaling and angiogenesis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

