

KD-Validated Anti-BCL2 Like 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1882

Product Overview

Name	KD-Validated Anti-BCL2 Like 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q07817
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human Bcl-XL
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 26 kDa, observed, 30 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q07817)

Gene Name	BCL2L1
About Gene Symbol	BCL2L1 (BCL2 Like 1) is a protein phosphatase involved in apoptosis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

