

KD-Validated Anti-Catalase Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1883

Product Overview



Name	KD-Validated Anti-Catalase Rabbit Monoclonal Antibody
Accession(Primary)	P04040
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Catalase
Species Reactivity	Human, Rat
Calculated MW	Predicted, 60 kDa, observed, 60 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P04040)

Gene Name	CAT
About Gene Symbol	CAT (Catalase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-catalase antibody (Cat#ABG1883). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-catalase antibody (Cat#ABG1883, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-catalase antibody (Cat#ABG1883). Catalase expression in wild type (WT) and catalase (CAT) shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-catalase antibody (Cat#ABG1883, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.