

KD-Validated Anti-E2F Transcription Factor 4 Rabbit

Rabbit monoclonal antibody

ABG1884

Product Overview

Name	KD-Validated Anti-E2F Transcription Factor 4 Rabbit Monoclonal Antibody
Accession(Primary)	Q16254
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human E2F4
Species Reactivity	Human
Calculated MW	Predicted, 44 kDa, observed, 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500 Flow Cytometry (FCM): 1:2,000

Target information (Q16254)

Gene Name	E2F4
About Gene Symbol	E2F4 (E2F Transcription Factor 4) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

