

KD-Validated Anti-NK2 Homeobox 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1889

Product Overview


Name	KD-Validated Anti-NK2 Homeobox 1 Rabbit Monoclonal Antibody
Accession(Primary)	P43699
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TTF1
Species Reactivity	Human
Calculated MW	Predicted, 39 kDa, observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P43699)


Gene Name	NKX2-1
About Gene Symbol	NKX2-1 (NK2 Homeobox 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genulin Biotechnologies LLC



Copyright ©2024 Genulin Biotechnologies LLC



Copyright ©2024 Genulin Biotechnologies LLC