

KD-Validated Anti-Ring Finger Protein 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1891

Product Overview

Name	KD-Validated Anti-Ring Finger Protein 2 Rabbit Monoclonal Antibody
Accession(Primary)	Q99496
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human RING2 / RING1B / RNF2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 38 kDa, observed, 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q99496)

Gene Name	RNF2
About Gene Symbol	RNF2 (Ring Finger Protein 2) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

