

KD-Validated Anti-Basigin (Ok Blood Group) Mouse

Mouse monoclonal antibody

ABG1899

Product Overview

Name	KD-Validated Anti-Basigin (Ok Blood Group) Mouse Monoclonal Antibody
Accession(Primary)	P35613
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Immunogen	Recombinant protein of human BSG
Species Reactivity	Human
Calculated MW	Predicted, 42 kDa, observed, 35-70 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P35613)

Gene Name	BSG
About Gene Symbol	BSG (Basigin (Ok Blood Group)) is a collagen extracellular matrix protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

