

KD-Validated Anti-Gephyrin Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1902

Product Overview

Name	KD-Validated Anti-Gephyrin Mouse Monoclonal Antibody
Accession(Primary)	Q9NQX3
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human GPHN
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 80 kDa, observed, 80 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9NQX3)

Gene Name	GPHN
About Gene Symbol	GPHN (Gephyrin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

